ANEROID BLOOD PRESSURE KIT COMFORT

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USER'S MANUAL

ANEROID BLOOD PRESSURE KIT COMFORT model A-40

CAUTION:

Please read all instructions before using this instrument to take your blood pressure. Follow all CAUTIONS to avoid damaging the instrument or hurting yourself. Talk to your doctor for information about your blood pressure.

TABLE OF CONTENTS


Unit Assembly	5
What Is Blood Pressure	6
How to Take a Reading	7
Quick Reference Guide	9
Care and Maintenance	10
Warranty	10
Technical Specifications	11

UNIT ASSEMBLY

Lay out the parts as pictured. The cuff should already be folded through the metal D-ring. If the cuff is not pre-folded, lay the cuff flat in front of you. Place the metal D-ring on your left. Pull the end of the cuff through the metal bar. If you fold the cuff correctly, the hook and pile material faces the outside of the cuff loop.

Look at the drawing to help you put the unit together.

1. Attach the gauge to the tubing on the left.
2. Attach the inflation bulb to the tubing on the right.
3. Attach the single end of the Y-tubing to the stethoscope chestpiece that is built-in into the cuff.
4. Attach the double end of the Y-tubing to the binaural. Attach the double end of the Y-tubing to the binaural.

NOTE: Model A-40 cuff fits arm circumference 22-36 cm.
For easy assembly follow diagram.

Connect each component securely to avoid injury and instrument damage.

WHAT IS BLOOD PRESSURE?


Blood pressure is the force exerted on the walls of your blood vessels as blood flows through them. The heart is like a pump. When it contracts, it sends a surge of blood through the blood vessels and pressure increases. This is called systolic pressure.

When your heart relaxes between beats, your blood pressure decreases. This is diastolic pressure.


When a doctor takes your blood pressure, he or she measures both your systolic and diastolic pressures and records them as numbers. For example, if your blood pressure reading is 126/76(126 over 76), your systolic is 126 and your diastolic is 76. The numbers are calculated in millimeters of mercury and recorded as 126/76 mmHg.

Blood pressure varies during the day. Factors influencing your blood pressure include physical activity, medications, and your emotional and physical condition. A single measurement does not provide a true picture of your blood pressure. You need to measure your blood pressure over a period of time.

IMPORTANT: CONSULT YOUR DOCTOR TO DETERMINE YOUR NORMAL BLOOD PRESSURE. ONLY A DOCTOR IS QUALIFIED TO DIAGNOSE AND TREAT HIGH BLOOD PRESSURE. PRACTICE THE MEASUREMENT PROCEDURES FOR THIS INSTRUMENT CAREFULLY. CHECK YOUR PROCEDURE WITH YOUR DOCTOR.

Heart Contracts
(Pressure Increases)


Heart Relaxes
(Pressure Decreases)

HOW TO TAKE A READING

1. Make sure you are seated with your feet flat and your back and arm supported. Rest your arm at heart level.


CAUTION: MAKE SURE ALL COMPONENTS ARE ASSEMBLED CORRECTLY AND SECURELY. FAILURE TO DO SO MAY RESULT IN INJURY OR AN INCORRECT READING.

2. Locate brachial artery by placing two fingers 2 cm above the bend of the elbow on the inside of the arm.


3. Put your left arm through the cuff loop. Place the bottom edge of the cuff approximately 1-2 cm above the elbow. Place the cuff so that the air-pipes were over the brachial artery.


4. The cuff should be wrapped so that it fits snugly and stays in place. You should be able to place one finger between the cuff and arm. Fasten the cuff by folding the end of the cuff over the D-ring and press the hook material firmly against the fuzzy pile material. Make sure the cuff sits correctly.


5. Insert stethoscope ear-pieces into ears. Make sure the chestpiece is in contact with your skin and is located over the brachial artery.

6. Rest your arm on a table. Relax your arm, and turn your palm upward. Ensure the cuff is same level as heart.


7. Place the gauge in front of you. Hold the inflation bulb in your right hand. Close the air release valve attached to the inflation bulb by turning it clockwise.


8. To inflate the cuff, rapidly squeeze the inflation bulb. If you KNOW your normal SYSTOLIC pressure, inflate the cuff until the gauge reaches approximately 30 mmHg ABOVE normal SYSTOLIC reading. If you are unsure or do not know your normal SYSTOLIC reading consult your doctor for the correct inflation level.

CAUTION: DO NOT INFLATE THE CUFF ABOVE 280 mmHg. YOU MAY INJURE YOURSELF OR DAMAGE THE INSTRUMENT.

9. Slowly open the air release valve by turning it counter-clockwise. Release the air at a rate of 5-6 mmHg per second.

10. Listen carefully for the appearance of sounds, watching the gauge needle. When you hear the first appearance of sound, it is your **SYS-TOLIC** blood pressure reading.

11. Continue to release the air at a rate of 5-6 mmHg per second. When you no longer hear any sounds, that is your **DIASTOLIC** blood pressure reading.

12. Turn the air release valve counter-clockwise to release the remaining air. Record your measurement in the record chart. You should record reading early, (not to forget values). Remove the stethoscope from your ears. Remove the cuff.

NOTE: If you want to take another reading, you can reinflate up to twice. Wait at least 5 minutes before measuring again, to avoid inaccurate readings due to the engorged blood vessels.

QUICK REFERENCE GUIDE

Before

1. **MEASURE** your blood pressure at the **SAME TIME** each day.
2. **RELAX** for at least 5 minutes before taking a measurement.
3. Remove tight fitting clothing from your upper arm.
4. **DO NOT** eat, smoke or exercise for at least 30 minutes before taking a measurement.

During

1. **PLACE CUFF** at the heart level.
2. **STETHOSCOPE CHESTPIECE** covers the **BRACHIAL ARTERY**.

3. INFLATE unit to proper level (30-40 mmHg above estimated systolic pressure).
4. DO NOT talk or move during a measurement.
5. DEFLATE unit at 5-6 mmHg per second.

After

1. WAIT 5 minutes before taking another measurement.

CARE AND MAINTENANCE

To protect your unit from damage, please AVOID washing or moistening the cuff, dropping the gauge, or hitting the surface of the stethoscope.

Your new blood pressure unit has been carefully checked to assure reliability and accuracy prior to shipment and use. However, as with any sensitive instrument subjected to repeated use, we recommend that your blood pressure gauge be checked periodically. When the cuff is fully deflated, the gauge needle must stay within the accuracy indicator zone. If the needle points outside of the accuracy indicator zone, the gauge will give inaccurate readings. In this case you have to bring the device to the nearest service centre for calibration (address mentioned in warranty card).

WARRANTY

The blood-pressure monitor is guaranteed for 2 years from the date of purchase. Warranty for the cuff and bulb is 1 year from the date of purchase. The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed by the seller (or by the service center representative) with the date of the purchase (repair), the check, instruction manual and safe packaging of the device.

- The warranty does not apply to wear out of cuff, rubber tubes and packaging box.
- The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.
- The device cannot be returned if it is without original packaging, with mechanical damages (scratches, stains and so on) or without any component.
- The device which was used and needs repair cannot be replaced with a new one.

In case there was violation of operating rules during warranty period, repair will be done at the expense of the consumer.

TECHNICAL SPECIFICATIONS OF MODEL A-40

Weight:	327 g
Size:	200 x 120 x 75 mm
Storage temperature:	-20°C to +70°C
Humidity:	10% - 90% relative humidity
Operation temperature:	5°C to 40°C
Measuring range:	20 to 300 mmHg
Measuring resolution:	2 mmHg
Accuracy:	± 3 mmHg
Air leakage:	$< 3 \pm 2$ mmHg/sec.

Accessories:

1. Sphygmomanometer.
2. Cuff (22-36 cm) with inlaid latex bladder.
3. Bulb and valve.
4. Stethoscope (the chestpiece is built-in into the cuff).
5. Soft bag.
6. Packaging.

Technical alterations reserved!

Whole service life without cuff - 10 years.

Whole service life of cuff - 2 years.

Checking

The device is calibrated at manufacturing checking. During usage and after repair it is recommended to check the device at least once a year in specialised labs.

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Wuxi Medical Instrument Factory Co., Ltd., No.43, Xixin Road, Zhangjing, Xibei Town,
Wuxi City, 214194 Jiangsu, China.  Lotus Global Co., LTD, 1 Four Seasons