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## DEAR GAMMA WRIST AUTOMATIC BLOOD PRESSURE MONITOR OWNER

Thank you for choosing Wrist Automatic Blood Pressure Monitor **Gamma**<sup>®</sup>TM model Active. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products Trademark **Gamma**<sup>®</sup>.

This system is a fully automatic, digital blood-pressure measuring device for use on the WRIST, it enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring. This monitor features Blood Pressure Level Indicator, Irregular Heartbeat Detection, as well as 30 memory cells.

Before starting to use this device please study the user's manual carefully. The user's manual offers all information you need to measure your blood pressure and pulse correctly. For all questions concerning the device please contact your local distributor or **Gamma**<sup>®</sup> service centre in your country.

### ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

Type BF applied part.

Read the instructions carefully before using this device.

# TABLE OF CONTENTS

---

<b>IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT</b> .....	<b>45</b>
What Is Blood Pressure? .....	45
Which Values Are Normal? .....	46
<b>ADVANTAGES OF WRIST AUTOMATIC BLOOD PRESSURE MONITOR MODEL ACTIVE</b> .....	<b>48</b>
Blood Pressure Level Indicator.....	48
Irregular Heartbeat Detection .....	49
Fuzzy Logic Measuring Technology .....	49
<b>GETTING READY FOR MEASUREMENT</b> .....	<b>50</b>
Safety Precautions .....	50
Description of the Blood Pressure Monitor .....	51
Insert or Replace Batteries .....	52
<b>MEASUREMENT PROCEDURE</b> .....	<b>53</b>
Fitting the Cuff.....	53
Taking a Single Measurement .....	54
<b>MEMORY FUNCTION</b> .....	<b>55</b>
Viewing the stored values .....	55
<b>ERROR MESSAGES/TROUBLESHOOTING</b> .....	<b>55</b>
<b>CARE AND MAINTENANCE</b> .....	<b>56</b>
<b>SYMBOL INFORMATION</b> .....	<b>57</b>
<b>TECHNICAL SPECIFICATIONS</b> .....	<b>58</b>
<b>WARRANTY</b> .....	<b>59</b>



# IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT

## What Is Blood Pressure?

Blood pressure is the pressure that your blood exerts on the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a «pump», sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.

There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteriae; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The

fluctuations are influenced by a number of factors - time of day, person's condition, physical or mental activity, environment, etc.

### **Day-Night Fluctuations of the Blood Pressure**

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic.

One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That's why the selfcontrol of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitalses (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

### **Which Values are Normal?**

---

Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age.

This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels, diabetes. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction.

Definitions and Classification of blood pressure levels according to 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension: Data in mmHg.

<b>Category</b>	<b>Systolic (mmHg)</b>		<b>Diastolic (mmHg)</b>
Optimal	< 120	and	< 80
Normal	120-129	and/or	80-84
High normal	130-139	and/or	85-89

Grade 1 hypertension	140-159	and/or	90-99
Grade 2 hypertension	160-179	and/or	100-109
Grade 3 hypertension	≥ 180	and/or	≥ 110
Isolated systolic hypertension	≥ 140	and	< 90

Isolated systolic hypertension should be graded (1, 2, 3) according to systolic blood pressure values in the ranges indicated, provided that diastolic values are < 90mmHg.

Source: The European Society of Hypertension and European Society of Cardiology Task Force Members. 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension. *J Hypertens* 2007; 25:1751-1762.

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

### ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental


exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor. When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

## ADVANTAGES OF WRIST AUTOMATIC BLOOD PRESSURE MONITOR MODEL ACTIVE

---

### Blood Pressure Level Indicator

---

Blood pressure level indicator is located along the left side in the display. The classification corresponds to ESH (European Society of Hypertension) ranges described in the table of the section «Which values are normal?». After the measurement there appears the bars in the left part of the display. This function helps you to self-orient in the measurement results.

## Irregular Heartbeat Detection

---

This function indicates allorhythmic heartbeating. If the symbol of IHD (♥) appears on the display that means that certain abnormality in heart beat frequency was detected during the measurement.

In this case, the result may deviate from your normal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol (♥) appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

### **Information for the doctor on frequent appearance of the Irregular Heartbeat Detector.**

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested. The symbol (♥) is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol (♥) appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice. The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

## Fuzzy Logic Measuring Technology

---

This unit uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure equivalent to the air pressure. This unit will determine the appropriate inflation level based on pressure oscillations, followed by cuff deflation.

During the deflation, the device will detect the amplitude and slope of the pressure oscillations and thereby determine for you the systolic blood pressure, diastolic blood pressure, and pulse.

# GETTING READY FOR MEASUREMENT

---


## Safety Precautions

---

1. Use this device ONLY for the intended use described in this manual.
2. Do NOT use accessories which are not specified by the manufacturer.
3. Do NOT use the device if it is not working properly or damaged.
4. Do NOT use under any circumstances on newborns.
5. This device does NOT serve as a cure for any symptoms or diseases. The data measured are for reference only. Always consult your doctor to have the results interpreted.
6. Keep the equipment and its flexible cord away from hot surfaces.
7. Do NOT apply the cuff to areas other than the place directed.
8. Proper maintenance and periodically calibration are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact local customer service for help.
9. This system is intended for individuals age 18 or above.
10. Do not use it for diagnosis of hypertension or testing on newborns, babies, young children or persons who cannot express their consent.
11. If you have been diagnosed with a severe arrhythmia or irregular heartbeat, atrial or ventricular premature beats or atrial fibrillation, measurements made with this instrument should only be evaluated after consultation with the doctor or your healthcare professionals.

**KEEP THIS MANUAL DURING THE WHOLE LIFETIME OF THE DEVICE**

## Description of the Blood Pressure Monitor

**Accessory:** Manual


Case

### SYMBOLS ON DISPLAY


WHO blood pressure classification

## Insert or Replace Batteries

---

1. Remove the battery cover.
2. Insert new batteries into the battery compartment as shown, taking care that the polarities (+) and (-) are correct.
3. Close the battery cover. Use only LR03, AAA batteries.

Disposal of old batteries to the authorized collecting party subject to the regulation of each individual territory.

### CAUTION

- Insert the batteries as shown in the battery compartment.
- When  (LOW BATTERY mark) blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.
-  (LOW BATTERY mark) does not appear when the batteries run out.
- Battery life varies with the ambient temperature and may be shorter at low temperatures.
- Remove the batteries if the device is not to be used for a long time.
- The batteries may leak and cause a malfunction.
- Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a shorter life.

## MEASUREMENT PROCEDURE

---

### ATTENTION

- Find time to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff.
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood pressure changes during the course of the day.
- Measurements should be done after a 5 minute rest to ensure accuracy.

### Fitting the Cuff

---

Select a pressure cuff size that corresponds to the circumference of your wrist (measured with a tight fit of the wrist). The size M (13,5 - 19,5 cm) is a size fitting the majority of people.

- a. Remove all eventual objects (e.g. wristwatch) from the wrist. Wrap the cuff over the wrist.
- b. The distance between the cuff and the hand should be as close as possible.
- c. Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight, no space should remain between the cuff and the wrist.
- d. Lay the wrist on a table, with the palm upwards. Support the wrist a little with a rest (cushion), so that the cuff rests at about the same height as the heart. Remain so for 2 minutes sitting quietly, before beginning with the measurement.

## Taking a Single Measurement

---

### ATTENTION

1. Sit down for at least 10 minutes before measuring.
2. Make sure the cuff is about the same height as the location of your heart.
3. Measure always on the same wrist (normally left).
4. Remain still and do not talk or move during the measurement.

**WARNING:** If the cuff is relatively lower than the heart, the obtained blood pressure value could be higher than the actual value. If the cuff is relatively higher than the heart, the obtained blood pressure value could be lower than the actual value. Always apply the pressure cuff before turning on the monitor.

1. Fasten the wrist cuff according to the instructions in “**ATTACHING THE WRIST CUFF.**”
2. Press the “ON/OFF” button. All icons appear two seconds on DISPLAY, then switch to measurement, and display “0” or last measurement record.

3. Start measurement, the cuff in the strap will automatically inflate. The mark will flash on LCD. Such measurements completed, LCD display measurement results.

## MEMORY FUNCTION

### Viewing the Stored Values

The device stores 30 most recent blood pressure test results.

Press the "ON/OFF" button 3 seconds, a memory read out the latest measurements, "ON/OFF" for the buttons (UP).

## TROUBLESHOOTING

If you have trouble in using the unit please check the following points first.

ERROR DISPLAY	POSSIBLE CAUSE	HOW TO CORRECT
Nothing is displayed. When you push the POWER button or  battery icon flash	No battery installation	Insert batteries
	Battery worn out	Replace with new batteries
	The polarities of batteries placed wrongly	Insert battery in the correct polarities
E1: Can't normally increase pressure	Check your wrist cuff for air leakage	Replace wrist cuff with new one
E3: Inflate pressure too high		Remeasure or send back to dealer to recalibrate
E2 E4: Hand shaking while measurement	Hand or body shaking while measurement	Keep static position and measure again

Battery icon on	Battery low power	Replace battery and measure again
The systolic pressure Value or diastolic Pressure value too high	1.The wrist cuff was held lower than your heart	Keep correct position and measure again
	2.The wrist cuff was not attached properly	
	3.You moved your body or spoke during measurement	
The systolic pressure value or diastolic pressure value too low	1.The wrist cuff was held higher than your heart	
	2.You moved your body or spoke during measurement	

The functionality of this device may be disturbed by the use of strong electromagnetic fields such as mobile phones or radio stations at close range, so we recommend that these devices be kept at least 1 m away. In cases where this is inevitable, please make sure that the device is working properly before use. If you are suffering from a heartbeat disorder (arrhythmia), evaluation of measurement results of the device can be given only after consultation with a physician.

## CARE AND MAINTENANCE

---

To avoid the monitor attracting dirt, dust or other contaminants, wash and dry your hands thoroughly before use.

### Cleaning

- To clean the monitor exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do NOT flush with water.
- Do NOT use organic solvents to clean the monitor.

- Do NOT wash the pressure cuff.
- Do NOT iron the pressure cuff.

### Disinfection

If required (for example, in a hospital setting), disinfection of the pressure cuff twice per week is recommended. Disinfect the interior of the pressure cuff (the surfaces that come into contact with the patient's skin) with the use of a soft cloth soaked in ethyl alcohol (75-90%) and wait for the pressure cuff to dry up. Cleaning of the pressure cuff after 200 uses is recommended.

### Monitor Storage

- Storage condition: -20°C to 55°C (-4°F to 131°F), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.

## SYMBOL INFORMATION

SYMBOL	REFERENT
	Consult instructions for use
	Manufacturer
	Serial number
	Caution, consult accompanying documents
	Dispose of in accordance with the requirements of your country

	Type BF Equipment
	Class II
 0413	CE mark

### Periodic calibration of device

The accuracy of measurement instruments must be verified from time to time. For this reason, verification of the indication of statistical pressure periodically (once a year) and after technical maintenance and repair is recommended. More information can be obtained from the service center.

## SPECIFICATIONS

Measuring Method	Oscillometric Measurement
Indication	Digital LCD display
Measuring Range:	Pressure: (30~280) mmHg Pulse: (40~199) Beat/min
Accuracy:	Static Pressure: $\pm 3$ mmHg Pulse: $\pm 5\%$
Memory:	30 Memories
Power supply:	2x1.5V Batteries (LR03 or AAA)
Operating condition:	+5°C~+40°C. 30% RH~80% RH Atmospheric pressure: 80 kPa~106 kPa
Storage condition:	-20°C~+55°C. 10% RH~93% RH Atmospheric pressure: 50 kPa~106 kPa

Dimensions:	Not more than 66(D) X 60(H) X 28(W) mm
Weight:	Not more than 0,150 kg, excluding batteries
Classification	Type BF
Wrist circumference	(13.5~19.5) cm

\* Specifications may be changed without notice in the event of improvement being made.

1. DSTU EN 1060-1:2015 "Non-Invasive Sphygmomanometer." General Requirements. Paragraph 9.2. DSTU EN 1060-3:2015 "Non-Invasive Sphygmomanometer." Part 3. Paragraph 9.2.

Type of protection against electric shock: INTERNALLY POWERED EQUIPMENT.

2. Degree of protection against electric shock: TYPE BF APPLIED PART.

3. Mode of operation: CONTINUOUS OPERATION.

4. Equipment not suitable for category AP&APG equipment use in presence.

#### STATEMENT

The system might not meet its performance specifications if used or stored above or below the following temperature and humidity ranges:

Operating conditions: +5°C~+40°C. 30% RH~80% RH

Storage conditions: -20°C~+55°C. 10% RH~93% RH

## WARRANTY

Your Wrist Automatic Blood Pressure Monitor Model Active is warranted for 5 years from date of purchase. Warranty for the cuff is 1 year from the date of purchase. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device. The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed.

Declaration of electromagnetic compatibility can be found at the link:

<http://medhouse-group.com.ua/distribjucija-i-servis>

ISO  
13485

CE  
0413

EC REP

Lotus NL B.V., Koningin Julianaplein 10, 1e  
Verd, 2596AA, The Hague, Netherlands.

Shenzhen Pango Electronic Co., Ltd No 25, 1st Industry Zone, Fenghuang Road, Xikeng Village, Henggang Town, Longgang District, Shenzhen City, Guangdong Province, China.

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