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DEAR GAMMA AUTOMATIC BLOOD PRESSURE MONITOR OWNER

Thank you for choosing Automatic Blood Pressure Monitor **Gamma**® TM model Plus. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Trademark **Gamma**®.

This system is a fully automatic, digital blood-pressure measuring device for use on the upper arm, it enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring. This monitor features Blood Pressure Level Indicator, Irregular Heartbeat Detection, as well as 90 memory cells.

Before starting to use this device please study the user's manual carefully. The user's manual offers all information you need to measure your blood pressure and pulse correctly. For all questions concerning the device please contact your local distributor or **Gamma**® service center in your country.

ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

Type BF applied part.


Read the instructions carefully before using this device.

TABLE OF CONTENTS

IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT	44
What Is Blood Pressure?.....	44
Which Values Are Normal?.....	45
ADVANTAGES OF AUTOMATIC BLOOD PRESSURE MONITOR	
MODEL Plus	47
Blood Pressure Level Indicator.....	47
Irregular Heartbeat Detection	47
Fuzzy Logic Measuring Technology	48
GETTING READY FOR MEASUREMENT	48
Safety Precautions	48
Description of the Blood Pressure Monitor.....	49
Inserting the Batteries.....	50
Using a Mains Adaptor	51
Setting the Date and Time.....	51
MEASUREMENT PROCEDURE	52
Fitting the Cuff	52
Taking a Single Measurement	53
MEMORY FUNCTION	55
Viewing the Stored Values.....	55
Delete Memories	55
ERROR MESSAGES/TROUBLESHOOTING.....	55
CARE AND MAINTENANCE	56
SYMBOL INFORMATION	57
TECHNICAL SPECIFICATIONS.....	58
WARRANTY	59

IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT

What Is Blood Pressure?

Blood pressure is the pressure that your blood exerts on the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a «pump», sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.


There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteriae; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors - time of day, person's condition, physical or mental activity, environment, etc.

Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic.

One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That's why the self control of the blood

pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitals (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

Which Values are Normal?

Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age.

This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels, diabetes. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction.

Definitions and Classification of blood pressure levels according to 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension: Data in mmHg.

Category	Systolic (mmHg)		Diastolic (mmHg)
Optimal	< 120	and	< 80
Normal*3	120-129	and/or	80-84
High normal	130-139	and/or	85-89

Grade 1 hypertension	140-159	and/or	90-99
Grade 2 hypertension	160-179	and/or	100-109
Grade 3 hypertension	≥ 180	and/or	≥ 110
Isolated systolic hypertension	≥ 140	and	< 90

Isolated systolic hypertension should be graded (1, 2, 3) according to systolic blood pressure values in the ranges indicated, provided that diastolic values are < 90mmHg.

Source: The European Society of Hypertension and European Society of Cardiology Task Force Members. 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension. *J Hypertens* 2007; 25:1751-1762.

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor.

When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

ADVANTAGES OF AUTOMATIC BLOOD PRESSURE MONITOR MODEL PLUS

Blood Pressure Level Indicator

Blood pressure level indicator is located along the left side in the display. The classification corresponds to ESH (European Society of Hypertension) ranges described in the table of the section «Which values are normal?». After the measurement there appears the bars in the left part of the display: green zone - optimal blood pressure, yellow - elevated, orange - too high, red - dangerously high. This function helps you to self-orient in the measurement results.

Irregular Heartbeat Detection

This function indicates allorhythmic heartbeating. If the symbol of IHD (♥) appears on the display that means that certain abnormality in heart beat frequency was detected during the measurement.

In this case, the result may deviate from your normal blood pressure - repeat the measurement. In most cases, this is no cause for concern. However, if the symbol (♥) appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Irregular Heartbeat Detector.

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested.

The symbol (♥) is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol (♥) appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice. The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Fuzzy Logic Measuring Technology

This unit uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure equivalent to the air pressure. This unit will determine the appropriate inflation level based on pressure oscillations, followed by cuff deflation.

During the deflation, the device will detect the amplitude and slope of the pressure oscillations and thereby determine for you the systolic blood pressure, diastolic blood pressure, and pulse.

GETTING READY FOR MEASUREMENT


Safety Precautions

1. Use this device ONLY for the intended use described in this manual.
2. Do NOT use accessories which are not specified by the manufacturer.
3. Do NOT use the device if it is not working properly or damaged.
4. Do NOT use the equipment where aerosol sprays are being used, or where oxygen is being administered.
5. Do NOT use under any circumstances on newborns.
6. This device does NOT serve as a cure for any symptoms or diseases. The data measured are for reference only. Always consult your doctor to have the results interpreted.
7. Keep the equipment and its flexible cord away from hot surfaces.
8. Do NOT apply the cuff to areas other than the place directed.
9. Proper maintenance and periodically calibration are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact local customer service for help.
10. This system is intended for individuals age 18 or above.
11. Do not use it for diagnosis of hypertension or testing on newborns, babies, young children or persons who cannot express their consent.
12. If you have been diagnosed with a severe arrhythmia or irregular

heartbeat, atrial or ventricular premature beats or atrial fibrillation, measurements made with this instrument should only be evaluated after consultation with the doctor or your healthcare professionals.


KEEP THIS MANUAL DURING THE WHOLE LIFETIME OF THE DEVICE

Description of the Blood Pressure Monitor

Cuff:

Plus offers a wide range cuff with D-ring, measuring size: 22~32cm, suitable for most of arm circumferences.

Accessory:

Manual


Adapter


Bag

Display screen:

Inserting the Batteries

To insert batteries:

1. Remove the cover of the battery compartment, which is located in the back side of the device (see the illustration).
2. Insert the batteries (4 AA, 1.5 V), following the indicated polarity. (+ -)

ATTENTION

- After the warning appears, device is blocked till the batteries have been replaced.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.
- Please use «AA» Long-Life or Alkaline 1.5V Batteries. Do not use rechargeable batteries.

Using a Mains Adaptor

It is possible to operate Automatic Blood Pressure Monitor Model Plus with a mains adaptor (output 6 V DC/ 500 mA).


1. Connect the micro plug to the socket at the device right side.
2. Plug the mains adaptor (output 6 V DC/ 500 mA) to a 100-240 V electrical power socket.

NOTE:

- Please use the adaptor of Trademark **Gamma**[®]. Any other adaptor being used, the service center is not considered responsible for the warranty servicing/repair of the monitor.
- No power is taken from the batteries while the mains adaptor is connected to the device.
- Please consult service center if you have questions relating to the mains adaptor.

Setting the Date and Time

1. Press and hold "SET" key until the year number displays and flashes on LCD to enter setting mode.
2. Press "MEM" key to adjust the year, then press "SET" key again to save your setting and enter the month setting mode.
3. Press "MEM" key to adjust the month. Follow the same steps to adjust date/hour/minute until setting completed.

NOTE: If the meter is idle for 1.5 minutes during the setting mode, it will switch off automatically.

MEASUREMENT PROCEDURE


ATTENTION

- Find time to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff.
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood pressure changes during the course of the day.
- Measurements should be done after a 5 minute rest to ensure accuracy.

Fitting the Cuff


Select a pressure cuff size that corresponds to the circumference of your shoulder (measured with a tight fit in the middle of the shoulder). The size M (XX-XX cm) is a size fitting the majority of people.

- a. Plug the cuff connector to the cuff socket on the device left side.

- b. Assemble the cuff as shown below. The smooth surface should be inside of the cuff loop and the metal D-ring should not touch your skin.
- c. Stretch your left (right) arm in front of you with your palm facing up. Slide and place the cuff onto your arm to let the air tube and artery mark region toward the lower arm.

Wrap and tighten the cuff above your elbow. The artery mark on the edge of the cuff should be approx. 0.8" to 1.2" (2 cm to 3 cm) above your elbow. Align the tube over the main arteries on the inside.

- d. Leave a little free space between the arm and the cuff: you should be able to fit two fingers between them. Remove all clothing covering or constricting the measurement arm.
- e. Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

Taking a Single Measurement

ATTENTION

1. Sit down for at least 10 minutes before measuring.
2. Place your elbow on a flat surface. Relax your hand with the palm facing up.
3. Make sure the cuff is about the same height as the location of your heart.
4. Remain still and do not talk or move during the measurement.

WARNING:

If the cuff is relatively lower than the heart, the obtained blood pressure value could be higher than the actual value. If the cuff is relatively higher than the heart, the obtained blood pressure value could be lower than the actual value.

Always apply the pressure cuff before turning on the monitor.

- a. Press ON/OFF button, pump begins to inflate the cuff. The rising pressure in the cuff is shown on the display.

b. Pressure in cuff increases to working pressure. Then pump stops and pressure falls.

NOTE: In case the inflation pressure is not sufficient to determine measurement result, the monitor automatically re-inflates to a higher level.

- c. The heart symbol (♥) starts to flash every heart beat once pulse is detected.
- d. After the measurement, systolic pressure, diastolic pressure, pulse rate and color bars classified from ESH (European Society of Hypertension) standards shows on the display.
- e. If irregular heartbeat was detected during the measurement, LCD displays the (♥) icon to remind users of heartbeat irregularity.

DISCONTINUING A MEASUREMENT


If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

MEMORY FUNCTION

Viewing the Stored Values

The device stores 90 most recent blood pressure test results along with respective dates and times in the monitor memory for each user.

When the monitor is turned off, press «MEM» button, a memory reading of the latest measurement will display. To see the next page (old measurements) press «SET».

Delete Memories

To delete memory, hold «MEM» button for five seconds, the LCD will display «no» - this means that all the memory has been

deleted. Troubleshooting

If you have trouble in using the unit please check the following points first.

ERROR DISPLAY	POSSIBLE CAUSE	HOW TO CORRECT
Nothing is displayed. When you push the POWER button or  battery icon flash	No battery installation	Insert batteries
	Battery worn out	Replace with new batteries
	The polarities of batteries placed wrongly	Insert battery in the correct polarities
E0	Testing pressure with Clothes	Take off the clothes and test it gain
	The marker of arm cuff  is not in alignment with the arm's artery.	Align the marker of arm cuff with the arm's artery, and measure again.
E1: Can't normally increase pressure	Check your cuff for air leakage	Replace the cuff with new one
E3: Inflate pressure too high		Remeasure or send back to dealer to recalibrate
E2 E4: Hand shaking while measurement	Hand or body shaking while measurement	Keep static position and measure again
Battery icon on	Battery low power	Replace battery and measure again
The systolic pressure Value or diastolic Pressure value too high	1.The cuff was held lower than your heart	Keep correct position and measure again
	2.The cuff was not attached properly	
	3.You moved your body or spoke during measurement	
The systolic pressure value or diastolic pressure value too low	1.The cuff was held higher than your heart	Keep correct position and measure again
	2.You moved your body or spoke during measurement	

The functionality of this device may be disturbed by the use of strong

electromagnetic fields such as mobile phones or radio stations at close range, so we recommend that these devices be kept at least 1 m away. In cases where this is inevitable, please make sure that the device is working properly before use. If you are suffering from a heartbeat disorder (arrhythmia), evaluation of measurement results of the device can be given only after consultation with a physician.

CARE AND MAINTENANCE

To avoid the monitor attracting dirt, dust or other contaminants, wash and dry your hands thoroughly before use.

Cleaning

- To clean the monitor exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do NOT flush with water.
- Do NOT use organic solvents to clean the monitor.
- Do NOT wash the pressure cuff.
- Do NOT iron the pressure cuff.

Disinfection

If required (for example, in a hospital setting), disinfection of the pressure cuff twice per week is recommended.

Disinfect the interior of the pressure cuff (the surfaces that come into contact with the patient's skin) with the use of a soft cloth soaked in ethyl alcohol (75-90%) and wait for the pressure cuff to dry up.

Cleaning of the pressure cuff after 200 uses is recommended.

Monitor Storage

- Storage condition: -20°C to 55°C (-4°F to 131°F), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.

SYMBOL INFORMATION

SYMBOL	REFERENT
	Consult instructions for use

	Manufacturer
	Serial number
	Caution, consult accompanying documents
	Dispose of in accordance with the requirements of your country
	Type BF Equipment
	Class II
	CE mark

Periodic calibration of device

The accuracy of measurement instruments must be verified from time to time. For this reason, verification of the indication of statistical pressure periodically (once a year) and after technical maintenance and repair is recommended. More information can be obtained from the service center.

SPECIFICATIONS

Measuring Method	Oscillometric Measurement
Indication	Digital LCD display
Measuring Range:	Pressure: (30~280) mmHg Pulse: (40~199) Beat/min
Accuracy:	Static Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Memory:	90 Memories

Cuff Size	22-32 cm
Power supply:	4x1.5V Batteries(LR6 or AA)
Operating condition:	+5°C~+40°C. 30% RH~80% RH Atmospheric pressure: 80 kPa~106 kPa
Storage condition:	-20°C~+55°C. 10% RH~93% RH Atmospheric pressure: 50 kPa~106 kPa
Dimensions:	Not more than 130(W) X 96(H) X 65(D) mm
Weight:	Not more than 0,380 kg, excluding batteries
Classification	Type BF

* Technical alterations are possible!

This device has been tested to meet the electrical and safety requirements of:

IEC/EN 60601-1,
IEC/EN 60601-1-2

Reference to Standards:

- DSTU EN 1060-1:2015 “Non-Invasive Sphygmomanometer.” General Requirements. Paragraph 9.2. DSTU EN 1060-3:2015 “Non-Invasive Sphygmomanometer.” Part 3. Paragraph 9.2.
- EN 1060-1/-3, NIBP-requirements
- IEC60601-1 General requirement for safety
- IEC60601-1-2 Requirements for EMC
- EN1060-4, NIBP clinical investigation
- AAMI/ANSI/IEC 80601-2-30, ANSI/AAMI/ISO 81060-2, NIBP requirements


WARRANTY

Your Automatic Blood Pressure Monitor Model Plus is warranted for 5

years from date of purchase. Warranty for the cuff and adapter is 1 year from the date of purchase. The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed. It needs to contain stamp of service center or selling point with date of purchase, user manual, undamaged original device package.

- The warranty does not apply to batteries and package.
- The warranty does not apply damage caused by improper handling, accidents, non-compliance with the operating instructions or attempts to open and/or repair the device.
- Devices without original package, with any mechanical damage (scratches, stains, etc.) and incomplete configuration can not be returned.
- If the device is used, and must be repaired, it can not be replaced with a new one.
- The repair is carried out at the expense of consumers in case, if device was used in violation of regulations during the warranty period.

Declaration of electromagnetic compatibility can be found at the link:
<http://medhouse-group.com.ua/distribjucija-i-servis>

B_PLUS_V3_121218

 Shenzhen Pango Electronic Co., Ltd No 25, 1st Industry Zone, Fenghuang Road,
Xikeng Village, Henggang Town, Longgang District, Shenzhen City, Guangdong Province, China.

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