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AUTOMATIC BLOOD PRESSURE MONITOR GAMMA™ MODEL SMART

INSTRUCTION BOOK

Thank you for choosing our Automatic Blood Pressure Monitor **Gamma®**™ designed on the basis of the most advanced technologies. The device you are holding now is your reliable assistant. You can measure your blood pressure with this monitor as accurately and precisely as an experienced doctor or a skilled nurse would do it.

We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of **Gamma®** Trade Mark.

Before starting to use this device please study the instruction book carefully. The instruction book offers all information you need to correctly measure your blood pressure and pulse.

Regular blood pressure control allows for reliable and easy health monitoring and for effective combating possible complications caused by blood pressure abnormalities.

ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

Type BF applied part.


Read the instructions carefully before using this device.

TABLE OF CONTENTS

IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT	46
What Is Blood Pressure?.....	46
Which Values Are Normal?.....	47
ADVANTAGES OF AUTOMATIC BLOOD PRESSURE MONITOR	
MODEL SMART	49
Blood Pressure Level Indicator.....	49
Irregular Heartbeat Detection	49
Fuzzy Logic Measuring Technology	50
GETTING READY FOR MEASUREMENT	50
Safety Precautions	50
Description of the Blood Pressure Monitor.....	51
Inserting the Batteries.....	52
Using a Mains Adaptor	53
Cuff Connectors.....	53
Setting the Date and Time	54
MEASUREMENT PROCEDURE	54
Fitting the Cuff	55
Measuring Procedure	56
Reading measurement results.....	56
Discontinuing a measurement.....	57
MEMORY FUNCTION	57
ERROR MESSAGES/TROUBLESHOOTING.....	58
CARE AND MAINTENANCE	59
SYMBOL INFORMATION	60
TECHNICAL SPECIFICATIONS	60
WARRANTY	61

IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT

What Is Blood Pressure?

Blood pressure is the pressure that your blood exerts on the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a «pump», sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.


There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteriae; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors – time of day, person's condition, physical or mental activity, environment, etc.

Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic.

One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That's why the self-control of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitals (encephalon, heart, kidneys, blood vessels) breaks down. Without special

treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

Which Values are Normal?

Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age.

This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels, diabetes. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction.

Definitions and Classification of blood pressure levels according to 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension: Data in mmHg.

Category	Systolic (mmHg)		Diastolic (mmHg)
Optimal	< 120	and	< 80
Normal*3	120-129	and/or	80-84
High normal	130-139	and/or	85-89
Grade 1 hypertension	140-159	and/or	90-99

Grade 2 hypertension	160-179	and/or	100-109
Grade 3 hypertension	≥ 180	and/or	≥ 110
Isolated systolic hypertension	≥ 140	and	< 90

Isolated systolic hypertension should be graded (1, 2, 3) according to systolic blood pressure values in the ranges indicated, provided that diastolic values are < 90mmHg.

Source: The European Society of Hypertension and European Society of Cardiology Task Force Members. 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension. *J Hypertens* 2007; 25:1751-1762.

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor.

When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

ADVANTAGES OF AUTOMATIC BLOOD PRESSURE MONITOR MODEL SMART

Blood Pressure Level Indicator

Blood pressure level indicator is located along the left side in the display. The classification corresponds to ESH (European Society of Hypertension) ranges described in the table of the section «Which values are normal?». After the measurement there appears the bars in the left part of the display: green zone - optimal blood pressure, yellow - elevated, orange - too high, red - dangerously high. This function helps you to self-orient in the measurement results.

Irregular Heartbeat Detection

This function indicates allorhythmic heartbeating. If the symbol of IHD (♥) appears on the display that means that certain abnormality in heart beat frequency was detected during the measurement.

In this case, the result may deviate from your normal blood pressure - repeat the measurement. In most cases, this is no cause for concern. However, if the symbol (♥) appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Irregular Heartbeat Detector.

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested.

The symbol (♥) is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol (♥) appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice. The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Fuzzy Logic Measuring Technology

This unit uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure equivalent to the air pressure. This unit will determine the appropriate inflation level based on pressure oscillations, followed by cuff deflation.

During the deflation, the device will detect the amplitude and slope of the pressure oscillations and thereby determine for you the systolic blood pressure, diastolic blood pressure, and pulse.

GETTING READY FOR MEASUREMENT

Safety Precautions


1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or damaged.
4. Do **NOT** use the equipment where aerosol sprays are being used, or where oxygen is being administered.
5. Do **NOT** use under any circumstances on newborns.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured are for reference only. Always consult your doctor to have the results interpreted.
7. Keep the equipment and its flexible cord away from hot surfaces.
8. Do **NOT** apply the cuff to areas other than the place directed.
9. Proper maintenance and periodically calibration are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact local customer service for help.
10. This system is intended for individuals age 18 or above.
11. Do not use it for diagnosis of hypertension or testing on newborns, babies, young children or persons who cannot express their consent.
12. If you have been diagnosed with a severe arrhythmia or irregular heartbeat, atrial or ventricular premature beats or atrial fibrillation, measurements made with this instrument should only be evaluated after consultation with the doctor or your healthcare professionals.

KEEP THIS MANUAL DURING THE WHOLE LIFETIME OF THE DEVICE

Components of Your Blood Pressure Monitor Kit

The blood pressure monitor kit includes Automatic Blood Pressure Monitor Model Smart, the cuff, 4 batteries (AA), bag for storing of devices, Instruction Book, Warranty Card, Kit Box, Bag, Adapter.

Description of the Blood Pressure Monitor

Cuff:

Smart offers a wide range rigid cuff Comfort measuring size: 22~42 cm, suitable for most of arm circumferences.

Accessory:

Manual


Adapter


Bag

SYMBOLS ON DISPLAY

Insert or Replace Batteries

1. Remove the battery cover.
2. Insert new batteries into the battery compartment as shown, taking care that the polarities (+) and (-) are correct.
3. Close the battery cover. Use only LR6, AA batteries.

CAUTION

- Insert the batteries as shown in the battery compartment.
- When  (LOW BATTERY mark) blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.
-  (LOW BATTERY mark) does not appear when the batteries run out.
- Battery life varies with the ambient temperature and may be shorter at low temperatures.
- Remove the batteries if the device is not to be used for a long time.

- The batteries may leak and cause a malfunction.
- Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a shorter life.

NOTE: To prolong the batteries' life the device switches off automatically if no button is pressed for 1 minute. Otherwise you can switch it off by pressing the START button.

Using an Mains Adaptor

It is possible to operate Automatic Blood Pressure Monitor Model SMART with a mains adaptor (output 6 V DC/ 500 mA).


1. Connect the micro plug to the socket at the device right side.
2. Plug the mains adaptor (output 6 V DC/ 500 mA) to a 100-240 V electrical power socket.

NOTE:

- Please use the adaptor of Trademark **Gamma**[®]. Any other adaptor being used, the service center is not considered responsible for the warranty servicing/repair of the monitor.
- No power is taken from the batteries while the mains adaptor is connected to the device.
- Please consult service center if you have questions relating to the mains adaptor.

Cuff connection


Insert the cuff connector into the opening provided on the side of the device as shown in the picture.


Setting the Date and Time

1. Under Power Off. Press and hold “SET” button until the year number displays and flashes on LCD to enter setting mode. And press  key to

adjust the year. Confirm your selection by pressing “SET” button.

2. Following the same steps to choose year, month, day, hour, minutes, units conversion (mmHg/kPa), press  to change the value.

ATTENTION: At installation of unit conversion be attentive. Choose units of measure mmHg.

MEASUREMENT PROCEDURE

ATTENTION

- Find time to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff.
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood

pressure changes during the course of the day.


- Measurements should be done after a 5 minute rest to ensure accuracy.

Fitting the Cuff

Select a pressure cuff size that corresponds to the circumference of your shoulder (measured with a tight fit in the middle of the shoulder). The size M-L (22-42 cm) is a size fitting the majority of people.

1. Plug the cuff connector to the cuff socket on the device left side.
2. Extend and place your left (right) arm in front of you with your palm up. Slide and place the cuff onto your arm to let the air tube and artery mark region toward the lower arm.

Wrap and tighten the cuff above your elbow. The artery mark on the edge of the cuff should be approx. 0.8" to 1.2" (2 cm to 3 cm) above your elbow. Align the tube over the main arteries on the inside.

3. Leave a little free space between the arm and the cuff: you should be able to fit two fingers between them. Remove all clothing covering or constricting the measurement arm.
4. Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

Measuring Procedure

1. Press  or  button, the display shows "0", Inflation icon ▲ flashes on display when it begins to inflate. During the inflation, please do not move.

2. Automatic deflation after measurement, and display the blood pressure, heart rate, and blood pressure indicator.

- Mistake, display "E", please refer the instruction about troubleshooting.

- When cuff remove during the measurement, please fix the cuff and test again.

- Stop the measurement in emergency, please press  or  power. Automatic deflation at the same time.

3. Power off. Press  or , to cut off the power.



Reading Measurement Results

When the measurement has been concluded, the measured blood pressure and pulse values are now displayed.

The appearance of this symbol (♥) signifies that an irregular heartbeat was detected. This indicator is only a caution.

It is important that you be relaxed, remain still and do not talk during measurements.

On the graphic indicator the element showing level of arterial pressure on a scale of World Health Organization is displayed

Note: We recommend contacting your physician if you see this indicator frequently.

To prolong the batteries' life the device switches off automatically if no button is pressed for 1 minute. Otherwise you can switch it off by pressing the  or  button.

Discontinuing a Measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the  or  button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

MEMORY FUNCTION

1. Under power off. Hold  or  for 3 seconds, get into relative memory mode, it reads the average of latest 3 measurements, press "SET" to find more memory.

2. Press  or  to cut off the power.

- If forget to turn it off. Automatically power off in one minute.
- More than 90 measurements, it will deletes the earliest one automatically
- The stored memory will not lose after replacing battery.

3. Delete memory (all the memory will be deleted) Showing memory, hold  or  for 5 seconds, display «» means all

the memories have been deleted.

TROUBLESHOOTING

If you have trouble in using the unit please check the following points first.

ERROR DISPLAY	POSSIBLE CAUSE	HOW TO CORRECT
Nothing is displayed. When you push the POWER button or  battery icon flash	No battery installation	Insert batteries
	Battery worn out	Replace with new batteries
	The polarities of batteries placed wrongly	Insert battery in the correct polarities
E0	Testing pressure with Clothes	Take off the clothes and test it gain
	The marker of arm cuff  is not in alignment with the arm's artery.	Align the marker of arm cuff with the arm's artery, and measure again.
E1: Can't normally increase pressure	Check your cuff for air leakage	Replace cuff with new one
E3: Inflate pressure too high		Remeasure or send back to dealer to recalibrate
E2 E4: Hand shaking while measurement	Hand or body shaking while measurement	Keep static position and measure again
Battery icon on	Battery low power	Replace battery and measure again
The systolic pressure Value or diastolic Pressure value too high	1.The cuff was held lower than your heart	Keep correct position and measure again
	2.The cuff was not attached properly	
	3.You moved your body or spoke during measurement	
The systolic pressure value or diastolic pressure value too low	1.The cuff was held higher than your heart	Keep correct position and measure again
	2.You moved your body or spoke during measurement	

The functionality of this device may be disturbed by the use of strong electromagnetic fields such as mobile phones or radio stations at close range, so we recommend that these devices be kept at least 1 m away. In cases where this is inevitable, please make sure that the device is working properly before use. If you are suffering from a heartbeat disorder (arrhythmia), evaluation of measurement results of the device can be given only after consultation with a physician.

CARE AND MAINTENANCE

To avoid the monitor attracting dirt, dust or other contaminants, wash and dry your hands thoroughly before use.

Cleaning

- To clean the monitor exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do NOT flush with water.
- Do NOT use organic solvents to clean the monitor.
- Do NOT wash the pressure cuff.
- Do NOT iron the pressure cuff.

Disinfection

If required (for example, in a hospital setting), disinfection of the pressure cuff twice per week is recommended.

Disinfect the interior of the pressure cuff (the surfaces that come into contact with the patient's skin) with the use of a soft cloth soaked in ethyl alcohol (75-90%) and wait for the pressure cuff to dry up.

Cleaning of the pressure cuff after 200 uses is recommended.

Monitor Storage

- Storage condition: -20°C to 55°C (-4°F to 131°F), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.

SYMBOL INFORMATION

SYMBOL	REFERENT
	Consult instructions for use

	Manufacturer
	Serial number
	Caution, consult accompanying documents
	Dispose of in accordance with the requirements of your country
	Type BF Equipment
	Class II
 0413	CE mark

Periodic calibration of device

The accuracy of measurement instruments must be verified from time to time. For this reason, verification of the indication of statistical pressure periodically (once a year) and after technical maintenance and repair is recommended. More information can be obtained from the service center.

TECHNICAL SPECIFICATIONS

Measuring Method	Oscillometric Measurement
Indication	Digital LCD display
Measuring Range:	Pressure: (30~280) mmHg Pulse: (40~199) Beat/min
Accuracy:	Static Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Memory:	Memory for 90 measurements for 2 users
Power supply:	4x1.5V Batteries (LR6 or AA)
Operating condition:	+5°C~+40°C. 30% RH~80% RH Atmospheric pressure: 80 kPa~106 kPa
Storage condition:	-20°C~+55°C. 10% RH~93% RH Atmospheric pressure: 50 kPa~106 kPa

Dimensions:	Approx: 140(W) X 110(H) X 48(D) mm
Weight:	Approx: 450 g, excluding batteries
Classification	Type BF
Cuff size	22-42 cm

* Technical alterations are possible!

This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2

Reference to Standards:

- DSTU EN 1060-1:2015 "Non-Invasive Sphygmomanometer." General Requirements. Paragraph 9.2. DSTU EN 1060-3:2015 "Non-Invasive Sphygmomanometer." Part 3. Paragraph 9.2.
- EN 1060-1/-3, NIBP-requirements
- IEC60601-1 General requirement for safety
- IEC60601-1-2 Requirements for EMC
- EN1060-4, NIBP clinical investigation
- AAMI/ANSI/IEC 80601-2-30, ANSI/AAMI/ISO 81060-2, NIBP requirements

WARRANTY

Your Automatic Blood Pressure Monitor Model Smart is warranted for 5 years from date of purchase. Warranty for the cuff and adapter is 1 year from the date of purchase. The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed. It needs to contain stamp of service center or selling point with date of purchase, user manual, undamaged original device package.

- The warranty does not apply to batteries and package.
- The warranty does not apply damage caused by improper handling, accidents, non-compliance with the operating instructions or attempts to open and/or repair the device.
- Devices without original package, with any mechanical damage (scratches, stains, etc.) and incomplete configuration can not be returned.
- If the device is used, and must be repaired, it can not be replaced with a new one.
- The repair is carried out at the expense of consumers in case, if device was used in violation of regulations during the warranty period

Declaration of electromagnetic compatibility can be found at the link:
<http://medhouse-group.com.ua/distribjucija-i-servis>

ISO
13485

CE
0413

 Shenzhen Pango Electronic Co., Ltd No 25, 1st Industry Zone, Fenghuang Road, Xikeng Village, Henggang Town, Longgang District, Shenzhen City, Guangdong Province, China.

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